

NEW!

MONARK TREADMILL FOR MEDICAL AND REHABILITATION

The Monark Medical Treadmill is designed for medical and rehabilitation purposes. Low entry speed, easy access and long handrails ensure the best usability.

The Monark Novo display together with the possibility to control the treadmill from a third party device (metabolic cart or ECG), enables easy and precise operating.



PANAMA COMMUNICATION

“USER-FRIENDLY, PRECISE AND ROBUST”

MAIN FEATURES

- Designed for medical and therapeutic purposes
- Low entry to facilitate easy access
- Full length handrails with extra support to simplify standing up from wheelchairs
- Low initial speed of 0.1 km/h (0.06 mph)
- Easy to connect to external devices
- Medical Device class IIa
- Monark Novo Control Unit

THE BASICS

- Speed range 0.1-20 km/h (0.06 – 12.4 mph)
- Speed incremental of 0.1 km/h
- Maximal user weight 250 kg
- Self-lubricating belt
- Protocols: Trackmaster and CSAFE
- BLE Smart and ANT+
- Training mode Novo Display: *Speed, inclination, distance, heart rate, time*



For more information, visit monarkexercise.se/treadmill

